**SureFood® GMO QUANT 35S Soya (2 x 50 Reakt.)**

Art. Nr. S2028

Mai 2017

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	Blockcycler/LightCycler® 480	Rotorcyclers
Initial Denaturation (HOLD) Cycles	5 min, 95°C 45	1 min, 95°C 45
Denaturation Annealing/Extension (CYCLE)	15 sec, 95°C 30 sec, 60°C	10 sec, 95°C 15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup	Detection: End of extension phase Reporter Dye: FAM Quencher Dye: TAMRA Passive Reference: none	LightCycler® Channel: 530/610 oder F1/F2 Acquisition mode: Single in extension phase Rotor-Gene Q Reporter Dye: FAM (Green)
Detaillierte Informationen zur Einstellung bestimmter real-time PCR Geräte stehen auf der CONGEN-Homepage zur Verfügung: http://www.congen.de/unternehmen/download		

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Description

The test is developed for the relative estimation of 35S CaMV promoter DNA amount in only Soya containing products. Therefore the kit contains two PCR systems, one specific for the 35S CaMV promoter sequence and the other one specific for Soya (the reference gene). The detection of the 35S CaMV promoter is according to German Food Law § 64 LFGB. The real-time PCR assay can be used with established real-time PCR instruments (Roche LightCycler®, Rotor-Gene Q, ABI PRISM, Eppendorf realplex, BioRad CFX96, Agilent MxSeries etc.). Because different GMO events contain a different number of 35S CaMV promoter copies systematically deviations may appear.

Limit of Detection

The 35S Soya PCR has a limit of detection of ≤ 5 DNA-copies. The assay limit of detection depends on sample matrix, processing grade, DNA-preparation and DNA-content.

The limit of quantitation depends on the concentration of the sample DNA used in the analysis. For example, if 50,000 target-sequence copies of Soya specific reference gene are present, the relative quantitation limit for 35S Soya DNA is 0.1 %.

DNA-preparation

For DNA-preparation of raw material the use of SureFood® PREP Basic and for highly processed food and feed the use of SureFood® PREP Advanced is recommended.

Kit components and storage

Kit Code	Reagenz	Menge	Deckelfarbe
1	Soya Reaction Mix	1100 µl	Yellow
2	35S Reaction Mix	1100 µl	Yellow
3	Taq Polymerase	11 µl	Red
4	Dilution Buffer	1300 µl	White
5	Standard DNA	50 µl	Dark Blue
6	Positive Control (1% RR-Soya)	100 µl	Light Blue

Store all reagents at -20°C and protected from light.

Additionally required equipment and materials

- real-time PCR instrument
- real-time PCR consumables (plates, tubes, foils, capillaries, caps)
- pipettes with filter tips
- unpowdered disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

Protocol

1. Preparation of the master-mix

Calculate the total number of reactions needed (samples, controls and standards).

Reactions needed for the Soya detection:

For each run: 5 reactions for the standard curve
3 reactions for controls (1x no-template and 2x Positive Control)

For each sample: at least 1 reaction with each sample DNA

Reactions needed for the 35S detection:

For each run: 5 reactions for the standard curve
3 reactions for controls (1x no-template and 2x Positive Control)

For each sample: at least 1 reaction with each sample DNA

May 2017

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use. The tube of the Taq Polymerase should be kept at -20°C and not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Soya Reaction Mix or 35S Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	20 µl	220.0 µl

Mix each master-mix well and centrifuge shortly before use.

2. Setup

	Blockcycler/LightCycler® 480	Rotorcyclers
Initial Denaturation (HOLD) Cycles	5 min, 95°C 45	1 min, 95°C 45
Denaturation Annealing/Extension (CYCLE)	15 sec, 95°C 30 sec, 60°C	10 sec, 95°C 15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup	Detection: End of extension phase Reporter Dye: FAM Quencher Dye: TAMRA Passive Reference: none	LightCycler® Channel: 530/610 or F1/F2 Acquisition mode: Single in extension phase Rotor-Gene Q Reporter Dye: FAM (Green)
Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage: http://www.congen.de/en/company/downloads		

3. Preparation of the standard DNA dilutions

Dilute the Standard DNA (**Code 5**) in 1:10 steps in Dilution Buffer (**Code 4**) in order to prepare different DNA concentrations for the standard curves of the Soya reference gene and the 35S detection gene. Prepare 5 dilutions of the supplied Standard DNA with the supplied Dilution Buffer. Prepare 5 reaction tubes (labeled S1 to S5) and add 45 µl Dilution Buffer (**Code 4**) each. The following procedure is recommended:

standard	dilution	copy number per µl	final copy number per reaction*
S1	45 µl Dilution Buffer + 5 µl Standard DNA	100,000 copies	500,000 copies
S2	45 µl Dilution Buffer + 5 µl DNA of S1	10,000 copies	50,000 copies
S3	45 µl Dilution Buffer + 5 µl DNA of S2	1000 copies	5000 copies
S4	45 µl Dilution Buffer + 5 µl DNA of S3	100 copies	500 copies
S5	45 µl Dilution Buffer + 5 µl DNA of S4	10 copies	50 copies

***Note:** 5 µl of standard DNA are used for each calibration point. The final copy number per reaction is to be entered in the analysis software of the real-time PCR detection system.

If the diluted DNA standards (S1 to S5) are not immediately used, store them at -20°C. The dilutions should be stable up to two months. Before use allow the reagents to thaw, mix them on a vortex and centrifuge carefully before opening and use.

4. Preparation of the PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells or capillaries.
- Close the tube of the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells or capillaries and close them.
- Pipette 5 µl of the Positive Control and the standard dilutions into the designated tubes/wells or capillaries and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates or capillaries into the PCR instrument and start the run according to the setup.

Interpretation of results

The calculation for both reactions (**Soya, 35S**) has to be made separately. Mark the standards, the controls and the samples for the specific system (**35S**) and make the evaluation according to the usual analysis program recommended by the real-time PCR instrument manufacturer. Repeat the same procedure for the Soya reference gene system. The value for the slope of the standard curve has to be between -3.1 and -3.6 and the correlation coefficient $R^2 > 0.98$. In case of different values for the standard curve, it should not be used for calculation.

By using the calculated copy numbers for **35S** and **Soya** the relative GMO content of the sample DNA and the Positive Control can be determined in the following way (example):

Sample 35S	1350 copies	Positive Control 35S	400 copies
Sample Soya	45,000 copies	Positive Control Soya	28,000 copies

Divide the copy number of specific system by the copy number of the reference gene system and multiply by 100 to obtain the percentage.

35S Soya content = $\frac{\text{35S copy number}}{\text{Soya copy number}} * 100$

sample DNA **35S Soya** content = $\frac{1350}{45,000} * 100 = 3 \%$ sample DNA **35S Soya** content = 3 %

For the given example the numbers lead to a **35S Soya** content of 3.0 % for the sample and 1.4 % for the Positive Control with the same calculation.

For a final calculation the use of a correction factor K for the correction of run-to-run fluctuations is necessary. The correction factor is the relation from the true percentage value of the Positive Control (1 % GMO content) and the measured GMO percentage of the Positive Control. The factor is calculated in the following way:

$K = \frac{\text{true GMO percentage of Positive Control}}{\text{measured GMO percentage of Positive Control}}$

K (example) = $\frac{1 \%}{1.4 \%} = 0.7$

The calculated value for the sample is multiplied with K to obtain a corrected value.

GMO percentage sample = measured GMO percentage of sample * K

sample (example) = $3.0 \% * 0.7 = 2.1 \%$

For that example the **35S Soya** content is 2.1 %.

Note: A positive 35S signal may result from DNA derived from CaMV (cauliflower mosaic virus, original organism of the 35S promoter) which sometimes is present for example in spices containing Brassicaceae species. With the virus specific detection system SureFood® GMO SCREEN CaMV it is possible to verify the detection of a genetic modification.

May 2017

Product Information

- Validation Report
- Microsoft Excel template of calculation
(Download: www.congen.de/en/company/downloads)

Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

Distribution and ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com