**CONGEN**

# **SureFood<sup>®</sup> PREP Advanced**

Art. No. S1053

50 extractions

## **User Manual**

Efficient DNA preparation from highly processed food and feed

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**Table of contents** 

1	General Information .....	10
1.1	Description .....	10
1.2	Kit components and storage .....	10
1.3	Additionally required equipment and materials .....	10
2	Protocol 1.....	11
2.1	Principle according to protocol 1 .....	11
2.2	Preparations.....	11
2.3	Extraction protocol.....	11
2.4	Product Information .....	12
2.5	Technical Support .....	12
2.6	Distribution and ordering .....	12
3	Protocol 2.....	13
3.1	Principle according to protocol 2 .....	13
3.2	Preparations.....	13
3.3	Extraction protocol.....	13
3.4	Product Information .....	15
3.5	Technical Support .....	15
3.6	Distribution and ordering .....	15
4	Safety information .....	16

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1x Proteinase K	<b>(Code K)</b>	1x Elution Buffer (10 ml)	<b>(Code E)</b>
1x PCR grade water (1,1 ml)	<b>(Code H)</b>	1x Elution Buffer X (10 ml)	<b>(Code X)</b>
1x Lysis Buffer (60 ml)	<b>(Code L)</b>	1x Spin Filter, grün (50 x)	<b>(Code F)</b>
2x Binding Buffer (10 ml)	<b>(Code B)</b>	2x Spin Filter, klar (50 x)	<b>(Code S)</b>
1x Pre-Wash Buffer (60 ml)*	<b>(Code P)</b>	4x Receiver Tubes 2,0 ml, klar (50 x)	<b>(Code R)</b>
1x Wash Buffer (60 ml)*	<b>(Code W)</b>	1x Receiver Tubes 1,5 ml, klar (50 x)	<b>(Code T)</b>

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# 1 General Information

## 1.1 Description

This kit is intended to be used for the extraction of animal and plant DNA (deoxyribonucleic acid). It can be used with two different protocols.

The DNA extraction protocol 1 is designed for the sensitive extraction of animal and plant DNA from food according to regulation (EU) 1169/2011.

The DNA extraction protocol 2 is designed for the extraction of plant DNA from highly processed food and feed. Furthermore it can be used for samples that produce PCR inhibiting DNA extracts when extracted with protocol 1. Potentially occurring PCR inhibiting substances are efficiently removed in additional binding and purification steps.

## 1.2 Kit components and storage

1x Proteinase K	<b>(Code K)</b>	1x Elution Buffer (10 ml)	<b>(Code E)</b>
1x PCR grade water (1.1 ml)	<b>(Code H)</b>	1x Elution Buffer X (10 ml)	<b>(Code X)</b>
1x Lysis Buffer (60 ml)	<b>(Code L)</b>	1x Spin Filter, green (50 x)	<b>(Code F)</b>
2x Binding Buffer (10 ml)	<b>(Code B)</b>	2x Spin Filter, clear (50 x)	<b>(Code S)</b>
1x Pre-Wash Buffer (60 ml)*	<b>(Code P)</b>	4x Receiver Tubes 2.0 ml, clear (50 x)	<b>(Code R)</b>
1x Wash Buffer (60 ml)*	<b>(Code W)</b>	1x Receiver Tubes 1.5 ml, clear (50 x)	<b>(Code T)</b>

\* After adding ethanol (purity  $\geq$  96 %, puriss; not supplied with the kit)

The components are divided into two boxes.

The Proteinase K should be stored at  $-20^{\circ}\text{C}$ . All other reagents of the kit should be stored dry and at room temperature ( $14\text{-}25^{\circ}\text{C}$ ).

Before every extraction make sure that all components have room temperature. If there are any precipitates within the provided solutions, solve them by warming carefully (up to  $30^{\circ}\text{C}$ ).

## 1.3 Additionally required equipment and materials

- suitable equipment for sample comminution and homogenization
- micro balance and spatula for weighing the samples
- reaction tubes free from DNA and DNase 1.5 ml; 2.0 ml
- waterproof pen and tags for labeling the reaction tubes
- unpowdered disposable gloves
- pipettes with filter tips
- Vortex mixer
- Thermomixer/ heating block (up to  $65^{\circ}\text{C}$ )
- micro centrifuge (up to 12,000 rpm)
- ethanol (purity  $\geq$  96 %, puriss) for preparation of Pre-Wash Buffer and Wash Buffer
- disposal bags or waste bin

A flow chart for visualization of the extraction protocol is available at the CONGEN homepage:

<http://www.congen.de/en/company/download>.

## 2 Protocol 1

The DNA extraction protocol 1 is designed for the sensitive extraction of animal and plant DNA from food according to regulation (EU) 1169/2011.

### 2.1 Principle according to protocol 1

1. Preparation of the basic material
2. Lysis at 65°C
3. Pre-filtration and setting of optimal binding conditions
4. Binding of the nucleic acids on a Spin Filter
5. Purification of the bound nucleic acids
6. Drying of the Spin Filter
7. Elution of nucleic acids from the Spin Filter

### 2.2 Preparations

#### General

- Add 1 ml PCR grade water (**Code H**) to the Proteinase K (**Code K**) and mix thoroughly.
- Add 30 ml ethanol to the Pre-Wash Buffer (**Code P**) and mix thoroughly.
- Add 42 ml ethanol to the Wash Buffer (**Code W**) and mix thoroughly.

#### Before each preparation

- Preheating the Elution Buffer (**Code E**) - Transfer the needed amount of Elution Buffer (**Code E**) under calculation of a reserve volume into a reaction tube (not provided with the kit) and equilibrate to 65°C (The Elution Buffer is necessary for step 7).

### 2.3 Extraction protocol

#### 1. Preparation of the basic material

Transfer 100 mg\* homogenized sample material into a 2.0 ml reaction tube (not provided with the kit).

\* The sample amount depends on the specific sample matrix. For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

#### 2. Lysis of the basic material

Add 580 µl of Lysis Buffer (**Code L**) and 20 µl of Proteinase K (**Code K**) to the reaction tube and mix it briefly. Incubate on a heating block under continuously shaking for 60 minutes at 65°C.

**Note:** For strongly swelling samples it can be necessary to add additional volume of Lysis Buffer during the lysis.

#### 3. Pre-filtration and setting of optimal binding conditions

Centrifuge the sample lysate for 1 min at 12,000 rpm.

Transfer the liquid supernatant into a new 1.5 ml reaction tube (not supplied with the kit) and centrifuge again for 1 min at 12,000 rpm.

Place a green Spin Filter (**Code F – green, for pre-filtration**) into a clear 2.0 ml Receiver Tube (**Code R**). Take 400 µl supernatant from the last centrifugation step and transfer the liquid directly onto the Spin Filter.

Centrifuge the Spin Filter with the Receiver Tube for 1 min at 12,000 rpm. After centrifugation discard the Spin Filter.

#### 4. Binding of the nucleic acids on a Spin Filter

Add 250 µl Binding Buffer (**Code B**) to the filtrate and mix.

Place a clear Spin Filter (**Code S – clear, for binding**) into a new, clear 2.0 ml Receiver Tube (**Code R**). Transfer the complete solution onto the Spin Filter. Incubate at room temperature for 1 min.

Centrifuge the Spin Filter with the Receiver Tube for 1 min at 12,000 rpm. After centrifugation discard the filtrate and place the Spin Filter into a new, clear 2.0 ml Receiver Tube (**Code R**).

#### 5. Purification of the bound nucleic acids

Add 550 µl Pre-Wash Buffer (**Code P**) to the Spin Filter and centrifuge at 1 min for 12,000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

Add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge at 1 min for 12,000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

Once more add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge at 1 min for 12,000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

#### 6. Drying of the Spin Filter

Remove the residual ethanol by final centrifugation for 2 min at 12,000 rpm.

#### 7. Elution of nucleic acids from the Spin Filter

Place the Spin Filter into a clear 1.5 ml Receiver Tube (**Code T**) and add 50 µl of the preheated (65°C) Elution Buffer (**Code E**) directly onto the Spin Filter. Incubate on a heating block for 3 min at 65°C (no shaking).

Centrifuge for 1 min at 10,000 rpm. After centrifugation discard the Spin Filter.

The eluted DNA is ready-to-use for the PCR. The DNA can be stored for up to 24 hours at +4°C. For a storage time of more than 24 hours it should be kept at -20°C.

#### 2.4 Product Information

- Material Safety Data Sheet
  - Validation Report
  - Flow chart
- (Download: <http://www.congen.de/en/company/download>)

#### 2.5 Technical Support

For further questions please contact your distributor or send an e-mail to [sales@r-biopharm.de](mailto:sales@r-biopharm.de).

#### 2.6 Distribution and ordering

R-Biopharm AG  
An der neuen Bergstrasse 17,  
64297 Darmstadt, Germany  
Phone: +49 (0) 61 51 - 81 02-0  
Fax: +49 (0) 61 51 - 81 02-20  
E-Mail: [orders@r-biopharm.de](mailto:orders@r-biopharm.de)  
[www.r-biopharm.com](http://www.r-biopharm.com)

## 3 Protocol 2

The DNA extraction protocol 2 is designed for the extraction of plant DNA from highly processed food and feed. Furthermore it can be used for samples that produce PCR inhibiting DNA extracts when extracted with protocol 1. Potentially occurring PCR inhibiting substances are efficiently removed in additional binding and purification steps.

### 3.1 Principle according to protocol 2

1. Preparation of the basic material
2. Lysis at 65°C
3. Pre-filtration and setting of optimal binding conditions
4. Binding of the nucleic acids on a Spin Filter
5. Purification of the bound nucleic acids
6. Drying of the Spin Filter
7. First Elution of nucleic acids from the Spin Filter
8. Repeated setting of optimal binding conditions
9. Second binding of the nucleic acids on a Spin Filter
10. Second purification of the bound nucleic acids
11. Drying of the Spin Filter
12. Elution of nucleic acids from the Spin Filter for analysis

### 3.2 Preparations

#### General

- Add 1 ml PCR grade water (**Code H**) to the Proteinase K (**Code K**) and mix thoroughly.
- Add 30 ml ethanol to the Pre-Wash Buffer (**Code P**) and mix thoroughly.
- Add 42 ml ethanol to the Wash Buffer (**Code W**) and mix thoroughly.

#### Before each preparation

- Preheating the Elution Buffer X (**Code X**) and Elution Buffer (**Code E**) - Transfer the needed amount under calculation of a reserve volume of Elution Buffer X (**Code X**) and Elution Buffer (**Code E**) into a reaction tube (not provided with the kit) and equilibrate to 65°C. The Elution Buffer X (**Code X**) is necessary for step 7 and Elution Buffer (**Code E**) is necessary for step 12.

### 3.3 Extraction protocol

#### 1. Preparation of the basic material

Transfer 150 mg\* homogenized sample material into a 2.0 ml reaction tube (not provided with the kit).

\* The sample amount depends on the specific sample matrix. For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

#### 2. Lysis of the basic material

Add 580 µl of Lysis Buffer (**Code L**) and 20 µl of Proteinase K (**Code K**) to the reaction tube and mix it briefly. Incubate on a heating block under continuously shaking for 30 minutes at 65°C.

**Note:** For strongly swelling samples it can be necessary to add additional volume of Lysis Buffer during the lysis.

### 3. Pre-filtration and setting of optimal binding conditions

Centrifuge the sample lysate for 1 min at 12,000 rpm.

Place a green Spin Filter (**Code F – green, for pre-filtration**) into a clear 2.0 ml Receiver Tube (**Code R**). Take 400 µl supernatant from the last centrifugation step and transfer the liquid directly onto the Spin Filter.

Centrifuge the Spin Filter with the Receiver Tube for 1 min at 12,000 rpm. After centrifugation discard the Spin Filter.

### 4. Binding of the nucleic acids on a Spin Filter

Add 250 µl Binding Buffer (**Code B**) to the filtrate and mix the sample.

Place a clear Spin Filter (**Code S – clear, for binding**) into a new, clear 2.0 ml Receiver Tube (**Code R**). Transfer the complete solution onto the Spin Filter. Incubate at room temperature for 1 min.

Centrifuge for 1 min at 12,000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

### 5. Purification of the bound nucleic acids

Add 550 µl Pre-Wash Buffer (**Code P**) to the Spin Filter and centrifuge at 1 min for 12,000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

Add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge at 1 min for 12,000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

### 6. Drying of the Spin Filter

Remove the residual ethanol by centrifugation for 2 min at 12,000 rpm.

### 7. Elution of nucleic acids from the Spin Filter

Place the Spin Filter into a new, clear 2.0 ml Receiver Tube (**Code R**) and add 200 µl of the preheated (65°C) Elution Buffer X (**Code X**) directly onto the Spin Filter. Incubate on a heating block for 3 min at 65°C (no shaking).

Centrifuge for 1 min at 10,000 rpm. After centrifugation discard the Spin Filter.

### 8. Repeated setting of optimal binding conditions

Add 200 µl of Lysis Buffer (**Code L**) and 200 µl of Binding Buffer (**Code B**) to the 200 µl filtrate from step 7 and mix thoroughly.

### 9. Second binding of nucleic acids on a Spin Filter

Place a new, clear Spin Filter (**Code S – clear, for binding**) into a new, clear 2.0 ml Receiver Tube (**Code R**). Transfer the complete solution from step 8 onto the Spin Filter. Incubate at room temperature for 1 min.

Centrifuge for 1 min at 12,000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

**10. Second purification of the bound nucleic acids**

Add 550 µl Pre-Wash Buffer (**Code P**) to the Spin Filter and centrifuge at 1 min for 12,000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

Add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge at 1 min for 12,000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

**11. Drying of the Spin Filter**

Remove the residual ethanol by final centrifugation for 2 min at 12,000 rpm.

**12. Elution of nucleic acids from the Spin Filter for analysis**

Place the Spin Filter into a new, clear 1.5 ml Receiver Tube (**Code T**) and add 50 µl of the preheated (65°C) Elution Buffer (**Code E**) directly onto the Spin Filter. Incubate on a heating block for 3 min at 65°C (no shaking).

Centrifuge for 1 min at 10,000 rpm. After centrifugation discard the Spin Filter.

**Note:** The DNA can also be eluted with a higher volume (100 - 200 µl) of Elution Buffer, depending on the expected yield of DNA.

The eluted DNA is ready-to-use for the PCR. The DNA can be stored for up to 24 hours at +4°C. For a storage time of more than 24 hours it should be kept at -20°C.

**3.4 Product Information**

- Material Safety Data Sheet
  - Validation Report
  - Flow chart
- (Download: <http://www.congen.de/en/company/download>)

**3.5 Technical Support**

For further questions please contact your distributor or send an e-mail to [sales@r-biopharm.de](mailto:sales@r-biopharm.de).

**3.6 Distribution and ordering**

R-Biopharm AG  
An der neuen Bergstrasse 17,  
64297 Darmstadt, Germany  
Phone: +49 (0) 61 51 - 81 02-0  
Fax: +49 (0) 61 51 - 81 02-20  
E-Mail: [orders@r-biopharm.de](mailto:orders@r-biopharm.de)  
[www.r-biopharm.com](http://www.r-biopharm.com)

## 4 Safety information

### Lysis Buffer

Warning  
H319 H412  
P280 P273 P305+P351+P338

### Proteinase K

Danger  
H315 H319 H334 H335  
P280 P305+P351+P338

### Binding Buffer

Danger  
H225 H319 H336  
P210 P233 P305+P351+P338

### Pre-Wash Buffer

Warning  
H302 H412 EUH032  
P280 P273 P305+P351+P338

<b>H225:</b>	Highly flammable liquid and vapour.
<b>H302:</b>	Harmful if swallowed.
<b>H315:</b>	Causes skin irritation.
<b>H319:</b>	Causes serious eye irritation.
<b>H334:</b>	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
<b>H335:</b>	May cause respiratory irritation.
<b>H336:</b>	May cause drowsiness or dizziness.
<b>H412:</b>	Harmful to aquatic life with long lasting effects.
<b>EUH032:</b>	Contact with acids liberates very toxic gas.
<b>P210:</b>	Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
<b>P233:</b>	Keep container tightly closed.
<b>P273:</b>	Avoid release to the environment.
<b>P280:</b>	Wear protective gloves/protective clothing/eye protection/face protection.
<b>P305+P351+P338:</b>	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

For further information we offer a Material Safety Data Sheet. Please contact your distributor or send an e-mail to [sales@r-biopharm.de](mailto:sales@r-biopharm.de).