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# VitaFast<sup>®</sup> Folic Acid

## Brief information

Easy to use microbiological microtiter plate test for the quantification of total folic acid (enriched and natural folic acid) in food, animal feed and pharmaceutical products. All required reagents and standard are contained in the test kit. The kit is sufficient for 96 determinations incl. standard. Evaluation requires a microtiter plate reader (610 - 630 nm, alternatively at 540 - 550 nm).

**The VitaFast<sup>®</sup> Folic Acid test has been validated and certified with the licence no.100903 as a Performance Tested Method by the AOAC Research Institute.**

## Sample preparation

Liquid samples (added folic acid):  
sterile filtration and dilution

Solid samples (added folic acid):  
sample homogenization, extraction, centrifugation  
and dilution

Liquid and solid samples (added and native folic acid):  
sample homogenization, enzymatic treatment,  
extraction, centrifugation and dilution

Time requirement: Test conduction .....approx. 60 min  
Evaluation.....2 min

Incubation: 44 - 48 h in the dark at 37 °C (98.6 °F)

Standard range: 0.16 - 1.28 µg / 100 g (ml)

Limit of quantitation: 0.16 µg / 100 g (ml)

Limit of detection: 0.018 µg / 100 g (ml)

Recovery: 90 - 105 %

Intra-assay CV for standards: < 10 %

Inter-assay CV for standards: < 10 %

## 1. Principle of the test

The VitaFast® Folic Acid microtiter plate test is a microbiological method for the quantitative determination of total folic acid (added and natural folic acid) in food, animal feed and in pharmaceutical products. The microbiological test system is in accordance with international norms.

Folic acid is extracted from the sample and the extract is diluted. The folic acid assay medium and the diluted extract are pipetted into the wells of a microtiter plate which is coated with *Lactobacillus rhamnosus*. The growth of *Lactobacillus rhamnosus* is dependent on the supply of folic acid. Following the addition of folic acid as a standard or as a compound of the sample, the bacteria grow until the vitamin is consumed. The incubation is done in the dark at 37 °C (98.6 °F) for 44 - 48 h.

The intensity of metabolism or growth of *Lactobacillus rhamnosus* in relation to the extracted folic acid is measured as turbidity and compared to a standard curve. The measurement is done using a microtiter plate reader at 610 - 630 nm (alternatively at 540 - 550 nm).

## 2. Reagents provided

Each test kit contains sufficient materials for 96 determinations incl. standards. Each test kit contains:

- 1 x microtiter plate with 96 wells, coated with *Lactobacillus rhamnosus*
- 3 x redist., sterile water (30 ml) for the preparation of assay medium, standards and dilution of extracted samples
- 3 x folic acid assay medium (solid)
- 3 x folic acid standard (solid)
- 3 x folic acid buffer (liquid)
- 3 x adhesive foil (1 complete and 2 halves of adhesive foils, sufficient for 3 test runs)
- 1 x additional holder for microtiter strips

### Remark:

No quality guarantee can be assumed after expiration of the shelf life.

### 3. Required reagents and instruments, not provided

- **sterile bench** (recommended for sterile working)
- **microtiter plate reader** 610 - 630 nm (540 - 550 nm)
- incubator with dark chamber, 37 °C (98.6 °F)
- water bath heatable to 95 °C (203 °F)
- pH meter
- centrifuge > 8000 x g
- sterile tips for graduated micropipette 20 - 200 µl; 100 - 1000 µl
- sterile graduated centrifuge vials with screw cap (15 and 50 ml) and sterile reaction vials 1.5 or 2.0 ml
- 500 ml screw glass jar and volumetric flasks (100 and 1000 ml), 100 ml beaker
- sterile filters polyethersulfon 0.2 µm with syringe
- redist. or deionized water for sample extraction
- NaOH 2 mol / l (8 g NaOH ad 100 ml redist. or deionized water)
- NaOH 1 mol / l and 0,1 mol / l
- phosphate buffer (0.05 mol / l, 0.1 % ascorbate; pH 7.2): solve 7.8 g sodium dihydrogen phosphate dihydrate and 1 g sodium ascorbate in 1 liter redist. or deionized water, adjust pH to 7.2; prepare the buffer solution daily fresh

#### Reagents for the determination of total folate

- pig pancreatin (e. g. Sigma P1750)
- chicken pancreatin (e. g. R-Biopharm P2002)

### 4. Warning and precautions for the user

- the assay medium can evoke irritations of mucosa, eyes and skin
- after running the test the strips used must be disposed of according to regulations (e.g. autoclaved)

### 5. Storage instructions

Store the kit / reagents at 2 - 8 °C (35.6 - 46.4 °F).

**Use the prepared reagents (standard, medium) directly and reject them after the assay.**

## 6. Sample preparation

For the determination of **added folic acid** in nutrient enriched solid samples, a hot water extraction is usually sufficient. Liquid samples can be used directly after sterile filtration and dilution with sterile water from the test kit. For the determination of the **total folic acid (native and added)** content, the sample has to be treated with enzyme.

Samples should be stored protected from light at 4 °C (39.2 °F). Standards and samples should be run in triplicate. **Unknown sample matrices** should be analyzed **with two dilutions** of the sample extract. Sample extracts have to be used within one day and should be stored in the dark until analyzing.

**Sample extraction** is carried out with 1 g (ml) homogenized sample in 40 ml redist. or deionized water or extraction solution. This equals a **sample extraction dilution factor of 40**. This factor is already included in the standard curve (see Quality Assurance Certificate). **For low vitamin concentrations a sample weight of up to 5 g (ml) can be used** (this has to be considered in the evaluation).

### Notes:

- Where an enzymatic sample preparation is followed, it is recommended that a reagent blank should be carried out to ensure that no folic acid present in reagents (enzymes) contribute to the final results.
- Spike recovery studies may also be performed in order to identify the presence of any inhibitory compounds, which may negatively affect the result. An aliquot of spiking solution containing a known amount of folic acid should be added to the sample prior to extracting. Spiking levels should generally be approximately 100% of the expected vitamin content of the sample, in order to obtain an accurate result. When low spike recoveries are obtained, this may indicate presence of an inhibitory compound, and the sample should be diluted further where possible to eliminate the compound. Spiking standards are available from R-Biopharm.

Only sterile sample extracts or sterile dilutions thereof should be pipetted onto the microtiter plate. Dilutions have to be prepared with sterile water from the test kit. **Therefore after the sample extraction sterile working conditions and sterile consumables are necessary.** A sterile filtration of the sample or the sample extract is **always** necessary for:

- samples like fruit juices and fitness drinks, which are not heated during sample extraction (except when the sample is heated 30 min at 95 °C (203 °F) in a water bath)
- samples containing herbs and spices as well as honey and tea
- vitamin mixes, premixes, tablets (highly enriched samples, see 6.3) (except when the sample is heated 30 min at 95 °C (203 °F) in a water bath)
- samples with low vitamin concentrations that are highly coloured (the filtration step eliminates the colouring)
- if filtration is not possible due to solid particles or due to cloudiness, centrifugation should be carried out before the sterile filtration step (greater than 8000 x g for 5 min)

**The sterile filtration of the sample is not necessary, if the sample extraction is carried out at 95 °C (203 °F) for 30 min. Nevertheless the extracted samples have to be diluted with sterile water from the test kit (the assay medium has to be always filtered).**

#### **Example for the dilution factors of the sample extract:**

Solid sample with a labeled concentration of 125 µg / 100 g

The dilution of the extract should be in the middle of the standard curve. Therefore, the expected concentration is divided by the concentration of standard 2.

#### **Calculation:**

125 µg / 0.32 µg = 390  
 → dilution factor 400 (1 : 400)

#### **Dilution steps:**

- a) 1 : 10 → 0.1 ml sample extraction solution + 0.9 ml sterile water from the test kit
- b) 1 : 10 → 0.1 ml from a) + 0.9 ml sterile water from the test kit
- c) 1 : 4 → 0.25 ml from b) + 0.75 ml sterile water from the test kit

### **6.1. Added folic acid in liquid samples (multivitamin juices, fitness drinks)**

Add 1 ml sample into a 50 ml sterile centrifuge vial and fill up exactly to 40 ml with redist. or deionized water. Thereafter shake, sterile filter (alternatively: heat the sample 30 min at 95 °C in a water bath, thereafter chill down quickly below 30 °C (86 °F)) and, depending on the concentration range, further dilute in 1.5 ml (or 2.0 ml) sterile reaction vials with sterile water from the test kit.

### **6.2. Added folic acid in fruit gums and candies**

Weigh about 15 - 20 g fruit gums or candies in a 50 ml sterile centrifuge vial, add about 40 ml redist. or deionized water and solve the sample at 95 °C (203 °F) in water bath. Chill down quickly to below 30 °C (86 °F). Transfer the extraction solution with redist. or deionized water quantitatively into a 100 ml volumetric flask and fill up to the mark with redist. or deionized water. Transfer the volume corresponding to 1 g of test material into a 50 ml centrifuge vial and fill up to exactly 40 ml with redist. or deionized water, shake, sterile filter (alternatively: heat the sample 30 min at 95 °C in a water bath, thereafter chill down quickly below 30 °C (86 °F)) and, depending on the concentration range, further dilute in 1.5 ml (or 2.0 ml) sterile reaction vials with sterile water from the test kit.

**Example:** sample weight 17 g fruit gums

calculation:  $100 \text{ ml} / 17 \text{ g sample weight} = 5.88 \text{ ml} / \text{g}$

1 g sample is contained in 5.88 ml. Transfer 5.88 ml into a 50 ml centrifuge vial and continue as described above.

### **6.3. Added folic acid in capsules, pills, vitamin mixes**

First the capsule or pill weight has to be determined (average of 5 - 10 capsules or pills). Then pills are homogenized in a mortar or mixer. Capsules are cut open and extracted.

### 6.3.1 Sample preparation with 1 g sample size and pre-extraction

Weigh 1 g of pills, vitamin mix, premix or cut open capsule into a 500 ml screw glass jar, add about 400 ml phosphate buffer (0.05 mol / l; 0.1 % ascorbate; pH 7.2, freshly prepared, see section 3) and shake well. Extract for 30 minutes at 95 °C (203 °F) in a water bath. During extraction the glass jar has to be shaken well at least five times. Chill down quickly to below 30 °C (86 °F). Transfer the extraction solution with redist. or deionized water quantitatively into a 1000 ml volumetric flask and fill up to the mark with redist. or deionized water. Transfer 1 ml into a 50 ml sterile centrifuge vial and fill up exactly to 40 ml with redist. or deionized water. Thereafter shake, sterile filter (alternatively: heat the sample 30 min at 95 °C in a water bath, thereafter chill down quickly below 30 °C (86 °F)) and, depending on the concentration range, further dilute in 1.5 ml (or 2.0 ml) sterile reaction vials with sterile water from the test kit.

#### **Attention:**

For the calculation of the result, the pre-dilution factor of 1000 has to be considered (1 g up to 1000 ml). The dilution step 1 ml up to 40 ml is already included in the standard curve.

### 6.3.2 Sample preparation with 0.2 g sample size

Weigh 0.2 g of pills, vitamin mix, premix or cut open capsule into a 50 ml sterile centrifuge vial, add about 30 ml phosphate buffer (0.05 M; 0.1 % ascorbate; pH 7.2, freshly prepared, see section 3), shake well and fill up to exactly 40 ml with phosphate buffer. Extract for 30 minutes at 95 °C (203 °F) in a water bath. During extraction the vial has to be shaken well at least five times. It is important to make sure that the centrifuge vials are tightly closed. Chill down quickly to below 30 °C (86 °F). Thereafter centrifuge (5 min > 8,000 x g) and, depending on the concentration range, further dilute the clear supernatant in 1.5 ml (or 2.0 ml) sterile reaction vials with sterile water from the test kit.

**Note:** For the evaluation the weight of sample has to be considered.

#### **6.4. Added folic acid in cereals, baby food, bread, flour, powder and meat products**

Weigh 1 g homogenized sample into a 50 ml sterile centrifuge vial, add about 30 ml phosphate buffer (0.05 M; 0.1 % ascorbate; pH 7.2, freshly prepared, see section 3), shake well and fill up to exactly 40 ml with phosphate buffer. Extract for 30 minutes at 95 °C (203 °F) in a water bath. During extraction the vial has to be shaken well at least five times. It is important to make sure that the centrifuge vials are tightly closed. Chill down quickly to below 30 °C (86 °F), centrifuge (5 min > 8,000 x g) and, depending on the concentration range, further dilute the clear supernatant in 1.5 ml (or 2.0 ml) sterile reaction vials with sterile water from the test kit.

#### **6.5. Total folic acid content (native and added folic acid) in milk products, cereals, baby food, meat and meat products**

To extract the bound, native folic acid as well, or to determine it in non fortified samples, the sample has to be extracted with enzyme. Procedures with additional treatments using enzymes are described in the literature. The following extraction procedure has proved itself:

Weigh exactly 1 g (ml) homogenized sample and 20 mg pig pancreatin (10 mg chicken pancreatin for samples like cereals, products containing cereals, vegetables, fruits, yeasts, yeast products or liver) into a 50 ml centrifuge vial, add about 30 ml phosphate buffer (0.05 mol / l; 0.1 % ascorbate; pH 7.2, freshly prepared, see section 3) and shake. Fill up to exactly 40 ml with phosphate buffer. Incubate for 2 h at 37 °C (98.6 °F) in the dark (shake at times). Thereafter, heat the extract for 30 minutes at 95 °C (203 °F) in a water bath. During extraction the vial has to be shaken well at least five times. It is important to make sure that the centrifuge vials are tightly closed. Chill down quickly to below 30 °C (86 °F), centrifuge (5 min > 8,000 x g) and, depending on the concentration range, further dilute the clear supernatant in 1.5 ml (or 2.0 ml) sterile reaction vials with sterile water from the test kit.

**Note:** For low folic acid contents or low sample dilutions (< 1 : 8) the blank of the enzymatic treatment solution has to be considered in the evaluation.

## 7. Test implementation

### 7.1. Test preparation

The **bottle with sterile water**: push the coloured lid up, pull off right up to the glass rim and turn entire lid to remove it.

**Folic acid standards** should be dissolved and diluted freshly. Each dilution is sufficient for three wells.

- open the folic acid standard bottle, place the lid down with the opening facing upwards
- add **x ml (x = see quality assurance certificate and label standard bottle)** sterile water (from the test kit) to the standard bottle. Close the standard bottle with the lid and dissolve the standard by shaking = **standard concentrate**
- take 6 sterile vials (1.5 - 2.0 ml) and prepare from the dissolved standard concentrate a standard curve according to the following scheme:

standard curve* in µg / 100 g (ml)	sterile water in µl		standard concentrate in µl		total volume in µl
blank: 0	900	+	0	=	900
standard 1: 0.16	900	+	100	=	1000
standard 2: 0.32	400	+	100	=	500
standard 3: 0.64	300	+	200	=	500
standard 4: 0.96	200	+	300	=	500
standard 5: 1.28	100	+	400	=	500

\*The sample extraction dilution of 1 : 40 is already included in the standard curve.

The **folic acid assay medium** is sufficient for 6 microtiter strips. Open the assay-medium bottle and remove the desiccant using tweezers (discard the desiccant).

- add 10 ml sterile water from the test kit to the folic acid assay medium bottle
- add additionally **1 ml of the folic acid buffer** from the test kit to the folic acid assay medium bottle (sufficient for one preparation)
- close the assay medium bottle carefully and shake well
- heat the bottle in a water bath to 95 °C (203 °F) for 5 min while shaking at least twice; always make sure that the bottle is tightly closed
- quickly chill down to room temperature below 30 °C (86 °F)
- filter the medium through a 0.2 µm filter into a sterile 15 ml centrifuge vial

## 7.2. Test procedure

Only sterile samples which are diluted with sterile water from the test kit should be pipetted onto the microtiter plate.

- remove the **required strips** of the microtiter plate and place them into the additional holder. Return the unused strips together with the desiccant to the foil bag and seal it well, store at 2 - 8 °C (35.6 - 46.4 °F)

**Pipette first the assay medium and then standard or diluted samples, as followed:**

- pipette 150 µl folic acid assay medium into the wells
- pipette 150 µl standard or diluted sample into the assigned wells (flush the pipette tip with standard or sample solution)
- cover the strips/cavities with adhesive foil: pull off the protective layer of the foil, place the foil flat onto the strips, smoothing it down by hand, press the foil firmly onto the strips  
important: make sure the cavities are sealed airtight by smoothing down the foil over the cavities, take special care with the wells around the edges
- incubate at **37 °C (98.6 °F)** in the dark for **44 - 48 h** in an incubator

## 7.3. Measurement

- press down the adhesive foil once more, place the microtiter plate upside down on a table and dissolve the microorganisms thoroughly by shaking the plate on the surface of the desk
- invert the plate to the regular position and remove the adhesive foil diagonally, starting from the upper right; **the foil is strongly adhesive, so pulling it off the microtiter plate must be done with great care: hold the strips firmly in the frame with the other hand while you pull off the foil diagonally from the top right to the back**
- destroy any bubbles on the surface of liquid in the wells (by means of a pipette tip or a needle)
- measure the turbidity with a microtiter plate reader at 610 - 630 nm (alternatively at 540 - 550 nm).

### **Note:**

- after 44 - 48 h of incubation, the microtiter plate can be stored for max. 48 h in the refrigerator, thereafter the turbidity should be measured
- to avoid any time losses due to weekends or bank holidays, the microtiter plate can be evaluated after 60 h. It is recommended to use a timer to turn off the incubator after 44 - 48 h

## 8. Evaluation

A **4-parameter** evaluation is recommended, e. g. the RIDA<sup>®</sup>SOFT Win from R-Biopharm.

**The test evaluation is correct on condition that**

- OD blank < OD standard 1
- OD standard 5 > 0.6 OD

**The sample dilution factor of 40 is already included in the standard curve (see Quality Assurance Certificate).** In the below formula merely the further dilution factor of the extract and a differing sample weight need to be taken into consideration.

$$\begin{array}{l} \text{Folic acid} \\ \text{(in } \mu\text{g / 100 g)} \\ \text{(in } \mu\text{g / 100 ml)} \end{array} = \frac{\text{conc. standard curve} \times \text{dilution factor}}{\text{sample weight in g (ml)}}$$

### Example:

Sample weight: 1 g  
Sample extraction dilution: **1 : 40 (must not be considered)**  
Dilution of the sample extraction: 1 : 400 (has to be considered)  
Measured concentration from the standard curve: 0.32  $\mu\text{g}$  folic acid / 100 g

$$0.32 \times 400 / 1 = 128 \mu\text{g folic acid} / 100 \text{ g}$$

### Evaluation for capsules, pills, vitamin mixes

Folic acid in  $\mu\text{g}$  / pill or capsule =

$$\frac{\text{conc. standard curve} \times \text{dilution factor} \times \text{pill or capsule weight in g}}{\text{sample weight in g} \times 100}$$

For technical assistance in the U.S. and Canada  
R-Biopharm Inc. USA, Michigan  
Telephone: (877) 789-3033  
Fax: (269) 789-3070

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

## 9. Literature

- EN 14131 : 2003 Microbiological Determination of Folate
- Official Methods of Analysis (AOAC) 960.46
- AOAC 992.05 Folic Acid in Infant Formula, Microbiological Methods
- AOAC 944.12 Folic Acid in Vitamin Preparations, Microbiological Methods
- Schweizer Lebensmittelbuch SLMB Kapitel 62 / 11.2.1
- Amtliche Sammlung Untersuchungsverfahren (ASU) L 00.00-87
- Vitamin-Bestimmungen, Verlag Chemie GmbH Weinheim/Bergstr.; Strohecker und Henning
- The Vitamins, 2. Auflage, Vol. VII. Academic Press, New York/London (1967): V. Herbert und J.R.Bertino
- AACC Approved Methods of Analysis (2001) First supplement to 10th Ed., Eagen Press, St. Paul, MN, Method 86 - 47
- Official Methods of Analysis (2005) AOAC INTERNATIONAL, Gaithersburg, MD, Method 2004.05 (in press)
- Koontz et al. (2005) Comparison of Total Folate Concentrations in Foods Determined by Microbiological Assay at Several Experienced U.S. Commercial Laboratories. Journal of AOAC INTERNATIONAL Vol. 88, No.3, 805 - 813
- Pandrangi et al. (2004) Optimization of microbiological assay of Folic Acid and determination of folate content in spinach. International Journal of Food Science and Technology 2004, 39, 525 - 532

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is a registered trademark of the ifp Institut für Produktqualität GmbH.