RIDA® ABSORBANCE 96

User manual

All rights reserved.

This User manual may be used only for its intended purpose. It may not be reproduced in whole or in part or be translated into another language without our express, prior written consent.

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© 2021 R-Biopharm AG, Darmstadt

Technical changes, deviations in illustrations and errors reserved.

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6.2 Cleaning.....	49
6.3 Technical Support	52
6.4 Repairs.....	52
7 Performance data and technical data.....	53
8 Warranty	54
EC Declaration of conformity/CE mark.....	55

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User manual

Dear customer,

We are delighted that you have chosen the RIDA® ABSORBANCE 96. To take full advantage of the instrument's performance and to enjoy your analysis instrument for many years, please read these instructions for use carefully before installation and commissioning and operate the RIDA® ABSORBANCE 96 in accordance with the instructions. The operating safety and function of the instrument can only be guaranteed if both general safety regulations and accident prevention regulations of the legislator as well as the safety instructions in this manual are regarded. We accept no liability for any damage resulting from improper use or incorrect operation.

Ensure that the manual is accessible at all times and are read and understood by all persons operating the RIDA® ABSORBANCE 96. Keep this manual in a safe place so that you can access it at any time if necessary.

Contents

1 General information.....	35
Liability for function and damage	36
Intended use	36
Warnings.....	36
Information	37
2 Instructions for using the RIDA®ABSORBANCE 96.....	38
2.1 Measurement method	38
2.2 About the RIDA®ABSORBANCE 96.....	38
2.3 Information on consumables	40
2.4 Power consumption.....	40
2.5 Expected life time.....	40
3 Preparing the product for use.....	41
3.1 Unpacking, storage and transport	41
3.2 Checking the supplied parts.....	42
3.3 Working environment	42
3.4 Setting up and installation	42
4 Operation	44
4.1 RIDASOFT® Win.NET	44
4.2 Switching on.....	44
4.3 Carrying out the experiment.....	44
4.4 Finishing the experiment and switching off	46
5 Information required by the user	47
5.1 Explanation of the signal lights.....	47
5.2 Notes on troubleshooting	47
6 Maintenance and cleaning	49
6.1 Maintenance	49
6.2 Cleaning.....	49
6.3 Technical Support.....	52

6.4 Repairs.....	52
7 Performance data and technical data.....	53
8 Warranty	54
Explanation of symbols (device nameplate).....	55
EC Declaration of conformity/CE mark.....	55

1 General information

Please carefully compare the contents of this delivery with the enclosed delivery note, the package leaflet or the invoice. We recommend that you keep a copy of this document together with the instructions, so that you can quickly access information about the date and scope of delivery in case of future queries, re-orders or service work.

Please be sure to remove all small parts from the packaging material.

Before using the RIDA® ABSORBANCE 96 for the first time, please read this user manual carefully.

We recommend that it always be kept with the RIDA® ABSORBANCE 96.

Explanation of symbols used

In this manual, particularly important remarks are highlighted as followed:

	Warning! Follow the instructions to avoid risk in personal injury.
	Warning! Follow the instructions to avoid damage to the instrument.
	Warning, electrical hazard!
	Biological hazard!
	Safety information! Follow the instructions to ensure optimal instrument performance.

Liability for function and damage

Liability for the function of the RIDA®ABSORBANCE 96 shall in any case pass to the owner or operator if the device is improperly maintained, repaired, or modified by persons, who do not belong to the authorized service personnel or if it is handled in a way that does not comply with its intended use.

The service and operation of the product must be in accordance with this manual. R-Biopharm AG shall not be liable for damages resulting from non-observance of the information above. Warranty and liability conditions of the terms of sale and delivery of R-Biopharm AG are not extended by the information above.

Intended use

The RIDA®ABSORBANCE 96 is a microtiter plate reader for the automated readout of RIDASCREEN® ELISA tests from R-Biopharm AG and ELISA sales products. RIDASCREEN® ELISA tests and R-Biopharm sales products are carried out by trained laboratory personnel according to the relevant test kit description under the stated conditions for the specific sample material.

Warnings

Warning!

An incorrect working environment can lead to a reduction in service life, damage to the device, or measurement errors. Please note the explanations and warnings in chapter 3.3.

Safety information!

To avoid measurement errors, it is essential to ensure that the instrument is connected properly and that the experiments are carried out correctly. Please note the explanations and warnings in chapter 4.

Warning!

Incorrect cleaning of the RIDA®ABSORBANCE 96 can reduce its service life and can cause damage to the instrument. Please note the explanations and warnings in chapter 6.2.

Warning, electrical hazard!

To avoid the risk of electric shock, do not expose the RIDA®ABSORBANCE 96 to rain or excessive moisture and do not touch the instrument or USB cable with wet hands.

Information

Safety information!

Follow all safety instructions on the device and in the attached documents.

Safety information!

Follow all general precautions that apply to electrical instruments.

Disposal of used electrical and electronic equipment

The symbol on the product or its packaging indicates that this product is not to be treated as normal household waste. It must be disposed at a collection point for the recycling of electrical and electronic equipment. By contributing to the correct disposal of this product, you protect the environment and the health of your fellow human beings. The environment and health are endangered by incorrect disposal. Material recycling helps to reduce the consumption of raw materials. For further information on how to recycle this product, please contact your local authority or municipal waste disposal companies.

Working with biological and harmful material

The RIDA®ABSORBANCE 96 is not to be used for the measurement of biohazardous substances.

Always notice the manufacturer's hazard information for the substances to be measured.

The RIDA®ABSORBANCE 96 does not produce any toxic or harmful gases or substances. During the measurement make sure that there are no toxic or harmful substances in the microtiter plate.

2 Instructions for using the RIDA® ABSORBANCE 96

2.1 Measurement method

The RIDA® ABSORBANCE 96 is designed to carry out fast and sensitive absorbance measurements. It measures the optical density (OD) of samples at defined wavelengths.

Absorbance/Transmission

Absorbance refers to the amount of light absorbed by a medium. Absorption reduces transmission. Transmission is the ratio of incident to transmitted light. Accordingly, the degree of transmission is calculated as follows:

$$T = (I/I_0)$$

where I is transmitted light and I_0 is incident light.

Optical Density

Optical density is a measure of the attenuation of light radiation after it has passed through a medium. Optical density is the logarithmic quantity that describes the reciprocal of the transmittance T :

$$OD = \log(I_0/I)$$

where I is transmitted light and I_0 is incident light.

Optical density is the absorbance of the sample plus other attenuating effects such as scattered light due to turbidity. It is therefore necessary to avoid scattered light in order to measure absorbance correctly.

If no other attenuating effects are present, then absorbance = optical density.

2.2 About the RIDA® ABSORBANCE 96

The RIDA® ABSORBANCE 96 is an absorbance measuring instrument, i.e. a measuring instrument for reading, recording and further processing of the absorbance values of a 96-well microtiter plate. The RIDA® ABSORBANCE 96 has 96 detection units, which allow measurement without a scanning mechanism.

In the upper part of the measuring instrument are 4 light sources with different wavelengths (540, 450, 405, 630 nm) and 4 interference filters. The four signal lights on the top of the reader show, which measuring channel is used for the current measurement.

The microtiter plate is inserted manually into the slot of the device. There is a lock at the bottom of the slot and an arrow for visual inspection that is visible if the microtiter plate is inserted completely.


Fig. 1: Graphical representation of the product from the front, top.

1: Signal lights, 2: Slot for microtiter plate, 3: Lock, 4: Arrow.

On the back of the RIDA® ABSORBANCE 96 is a USB plug for connecting the instrument to a computer using the included Micro-USB cable. Furthermore, there is an eject button for pushing the microtiter plate out of the slot.


Fig. 2: Graphical representation of the RIDA® ABSORBANCE 96 from behind, top.
5: Eject button, 6: Micro-USB plug

2.3 Information on consumables

With the RIDA® ABSORBANCE 96 ANSI/SBS Standard 96 well microtiter plates can be readout.

When using microtiter plates, always note the specifications of the manufacturer of the microtiter plates. Not all microtiter plates of a manufacturer are the same in design, materials, or configuration. The temperature stability within the microtiter plate can depend on the type of microtiter plate used.

2.4 Power consumption

The RIDA® ABSORBANCE 96 is an analytical instrument with very low power consumption. It is operated by a Micro-USB cable and the total power consumption is generally less than 2.5 watts.

2.5 Expected life time

With a daily use of the RIDA® ABSORBANCE 96 of 4 hours the expected life time is 10 years.

3 Preparing the product for use

3.1 Unpacking, storage and transport

Unpacking

Remove the packaging material and carefully place the unit on a firm and level surface.

Check the parts supplied (see chapter 3.2).

Safety information!

Check the product for external damages.

Safety information!

Make sure that the rubber feet are completely attached to the underside of the device.

The RIDA®ABSORBANCE 96 is packed in a specially designed cardboard box. Keep the packaging material. If you return the device for repair, the original packaging material and the box must be used.

Storage and protection during the intervals of normal use

Warning!

Protect the RIDA®ABSORBANCE 96 from moisture and dust during prolonged storage.

Warning!

Do not store the RIDA®ABSORBANCE 96 next to heat sources and protect it from direct sunlight.

The storage temperature should be between 5 °C and 40 °C.

Transportation

Before transportation, make sure that the RIDA®ABSORBANCE 96 is disconnected from the computer and that it is turned off. Remove the cable from the instrument and make sure that there is no microtiter plate in the instrument. Depending on the transport distance, use the original packaging material. Make sure that the new location meets the requirements described in chapter 3.3.

Warning!

Always transport the product well packed to avoid possible transport damage.

Warning!

Before transport, make sure that there is no microtiter plate in the instrument.

3.2 Checking the supplied parts

When unpacking the product, please check that the following components are present:

- RIDA® ABSORBANCE 96
- Micro-USB cable
- Hex key
- Microfiber cloth

3.3 Working environment

Due to its small size, the RIDA® ABSORBANCE 96 can be operated very well in various working environments.

However, the following warnings must be regarded. Ignoring them may lead to a reduction in expected life time, damage to the device, or measurement errors:

Warning!

The RIDA® ABSORBANCE 96 should stand on a level surface, free from dust and vibrations.

Warning!

Do not operate the RIDA® ABSORBANCE 96 near heat sources or under direct sunlight.

Warning!

The ambient temperature should be between 5 °C and 40 °C.

Warning!

The maximum humidity should not exceed 90 %.

3.4 Setting up and installation

Step 1: Set up and make the RIDA® ABSORBANCE 96 ready for use

Warning!

Place the RIDA® ABSORBANCE 96 in a horizontal position and secured against falling.

Step 2: Connecting with the computer

Connect the RIDA®ABSORBANCE 96 to a computer using the Micro-USB cable. Connect the instrument directly to your computer and do not use an external USB hub. If the connection is successful, the four signal lights light up one after the other.

Warning!

Always use the included Micro-USB cable for the connection to your computer.

Warning!

The RIDA®ABSORBANCE 96 may only be operated at a USB port USB 2/ USB 3 with 5 VDC and a maximum of 3 A of a certificated computer (with certification mark of an approved testing laboratory). Use an adapter cable on USB 3.1 with type C plug connection, which ensures the profile 1: 5 V@2.0 A.

From Windows 98 drivers are available and have to be installed.

From Windows 8 drivers will be installed automatically (with internet connection).

If no connection between the RIDA®ABSORBANCE 96 and the computer is established, install the driver for the USB connection. You can find it at:

<https://www.ftdichip.com/Drivers/VCP.htm>

Step 3: Settings in the software

Open the RIDASOFT® Win.NET software and select "RIDA®ABSORBANCE 96" in the reader settings (File – Reader Setup). Make sure that the following settings have been made:

- Serial port: Serial port of the reader
- Baud rate: 115200
- Databits: 8
- Parity: Non
- Stop bits: 1

4 Operation

Before commissioning the RIDA® ABSORBANCE 96, you should carefully read and understand the entire manual to be familiar with the system.

4.1 RIDASOFT® Win.NET

The control of the RIDA® ABSORBANCE 96, presetting, data output, and analysis are controlled via a separate software solution (RIDASOFT® Win.NET).

Therefore, the software must be installed in order to use the instrument. To do this, follow the instructions in the Software user manual.

Start the RIDASOFT® Win.NET software.

4.2 Switching on

The RIDA® ABSORBANCE 96 requires a connection to the computer for power supply, to start the experiment, and for data transfer. Insert the Micro-USB cable into the RIDA® ABSORBANCE 96 and the USB port of the computer.

The RIDA® ABSORBANCE 96 starts automatically when connected to a computer via the Micro-USB cable and the RIDASOFT® Win.NET automatically establishes a connection.

When the RIDA® ABSORBANCE 96 is switched on, the signal lights on the top of the instrument light up. Each time the instrument is switched on an internal self-test is performed to ensure that there are no malfunctions.

4.3 Carrying out the experiment

For details on how to carry out experiments, please read the manual of the RIDASOFT® Win.NET. A brief description of the procedure is given below.

Before starting an experiment, make sure that the RIDA® ABSORBANCE 96 is properly connected to the computer and the RIDASOFT® Win.NET software.

The experiment consists of two steps:

1. Calibration

Start in the RIDASOFT® Win.NET software a new assay and choose, if you want to use a reference channel. Afterwards, a calibration in the measurement channels will be executed automatically.

This is indicated by a short illumination of the respective signal lights.

Warning!

Make sure that there is no microtiter plate in the slot during calibration. Only insert the microtiter plate into the slot of the RIDA® ABSORBANCE 96 **after successful calibration**.

2. Measurement

After successful calibration, you are requested to insert a microtiter plate. Now, insert the microtiter plate into the slot.

Pay attention to the alignment of the microtiter plate.

For orientation, a pictogram on the top of the RIDA® ABSORBANCE 96 indicates the correct alignment of the microtiter plate. (The pictogram is a schematic representation. The actual shape of the microtiter plate, especially the position of the rounded corner, may vary and differ from the illustration).

After insertion confirm in the RIDASOFT® Win.NET software that the microtiter plate has been inserted. Now, the measurement is carried out, which can be seen by a short illumination of the respective signal lights. Afterwards, the results are displayed in the RIDASOFT® Win.NET software.


Fig. 3: Insert the microtiter plate.

In order to guarantee a correct measurement result, please pay attention to the following warnings:

Safety information!

When inserting the microtiter plate, pay attention to its alignment.

Safety information!

Make sure that the microtiter plate was inserted in the slot completely. In that case the arrow at the bottom of the slot is visible completely.

**Warning!**

The bottom of the microtiter plate must be dry before it is inserted into the instrument. If the microtiter plate is wet at the bottom, dry it with a paper towel or similar.

**Safety information!**

Shocks of the instrument and of the microtiter plate can cause an unwanted signal. Therefore, neither the RIDA® ABSORBANCE 96 nor the microtiter plate should be touched during measurement.

**Safety information!**

Note that the system depends on a clean and dust-free environment. Therefore, do not wear powdered safety gloves.

**Safety information!**

Ensure that no direct sunlight or other strong light sources shine on the RIDA® ABSORBANCE 96 and especially into its slot.

4.4 Finishing the experiment and switching off

After completion of the experiment, the microtiter plate can be removed from the RIDA® ABSORBANCE 96. To do this, push the eject button on the rear side.


Fig. 4: Remove the microtiter plate.

Disconnect the Micro-USB cable to turn off the RIDA® ABSORBANCE 96.

5 Information required by the user

5.1 Explanation of the signal lights

Tab.1

Activity of the signal lights	Meaning
All signal lights flash simultaneously	Error, software shows error message
No signal light lights up	<ul style="list-style-type: none">• Device is off• Device is waiting for input• Device is defective
All signal lights light up briefly one after the other	Self-test after switching on the device
Signal light of the respective measuring channel flashes/illuminates continuously	Calibration/measurement in progress

5.2 Notes on troubleshooting

If there is an error, all four signal lights flash simultaneously. In this case, check the software to see which error is present.

In the event of the following errors, a new measurement can be carried out immediately after the error has been corrected:

Tab. 2

Type of error	Causes and correction of errors
Calibration failed	During calibration, there was a microtiter plate in the slot: Check if there is a microtiter plate in the slot and remove it before the instrument performs calibration.
	The inside of the slot is dirty: Clean the inside of the slot according to the cleaning instructions in chapter 6.2.
Interference due to ambient light, that is too bright	Too much ambient light has entered the slot: Make sure that the instrument is not exposed to too much ambient light, e.g. direct sunlight, during the measurement process.
Temperature out of specification	The ambient temperature was too high or too low: Make sure that the ambient temperature is between 5 °C and 40 °C.

In case of the following error, a new measurement cannot be performed directly. A restart of the RIDA® ABSORBANCE 96 is necessary first (remove the cable and reconnect it):

Tab. 3

Typ of error	Causes and correction of errors
USB power supply defective (< 450 mA or < 4 V)	Defective USB port: Use a different USB port on the computer or use a different computer.
	Using an USB hub: Connect the instrument directly to the computer, i.e. without an USB hub. If this is not possible use another USB hub.
	Defective Micro-USB cable: Contact your local distributor or R-Biopharm AG, Darmstadt.

If the following error occurs, the user cannot correct the error:

Tab.4

Type of error	Causes and correction of error
Irrecoverable hardware error	Irrecoverable hardware error: Contact your local distributor or R-Biopharm AG, Darmstadt.

6 Maintenance and cleaning

6.1 Maintenance

The RIDA®ABSORBANCE 96 is maintenance-free. Each time the instrument is switched on, an internal self-test is carried out to ensure that there are no malfunctions.

To check the linearity and accuracy of the measuring instrument for routine quality controls, a reference plate, e.g. from Hellma GmbH (item code 666-R113) can be used.

There are no parts within the RIDA®ABSORBANCE 96 that can be serviced by the customer. It is only necessary to ensure that the device is kept clean.

6.2 Cleaning

Remove the Micro-USB cable from the RIDA®ABSORBANCE 96 before cleaning.

Warning!

Make sure before cleaning that there is no microtiter plate in the RIDA®ABSORBANCE 96.

Cleaning the housing

The surfaces of the housing should be cleaned regularly. For this purpose, you can use e.g. a cloth or a sponge lightly wetted in water. For heavier soiling, clean the surface of the housing with a mild soap solution diluted with water or glass cleaner and then wipe with a lightly moistened cloth or sponge to remove any residue.

Do not use scouring agents.

Cleaning the slot

To clean the slot of the RIDA®ABSORBANCE 96, first open the four screws on the bottom of the instrument (see Figure 5) using the hex key delivered separately.


Fig. 5: Bottom side of the RIDA® ABSORBANCE 96.

Warning!

Do not loosen any screws that are not explicitly marked in the manual, as this can lead to a malfunction and may void the instrument's warranty.


Fig 6: Pull the RIDA® ABSORBANCE 96 apart into two parts.

Carefully pull apart the top and the bottom of the RIDA®ABSORBANCE 96 to unplug the plug connection on the back of the instrument.

Warning!

To avoid damaging the rigid plug connection between the top and bottom of the RIDA®ABSORBANCE 96, grasp the instrument at the back during the pull movement.

Avoid any tilting movement.

Use a dry microfiber cloth or a bellows to carefully remove dust or dirt from the surface. If this is not enough, the microfiber cloth can be slightly soaked with 80% ethanol to carefully remove dust or contamination from the surface.

Do not allow water or other liquids to get inside the RIDA®ABSORBANCE 96. If this happens, return the device to the manufacturer for inspection.

Warnings for cleaning

When cleaning, pay attention to the following warnings. The non-observance of these warnings may result in damage to the user or to the instrument. This can lead to a reduction of the service life or in measurement errors:

Warning!

Always avoid spraying liquid directly onto the surfaces of the RIDA®ABSORBANCE 96. This is especially important for the inside of the slot, where there are optical elements that are very sensitive and essential for the functionality of the instrument.

Warning!

Never clean the inside of the slot with sharp, abrasive and solvent-containing agents.

Biological hazard!

Always wear gloves during cleaning operations that may be associated with contact with biological or generally hazardous materials or liquids.

6.3 Technical Support

In case of problems with the RIDA® ABSORBANCE 96, please contact the technical service or the local responsible person.

Biological hazard!

It is your responsibility to decontaminate the instrument and all accessories before servicing and before returning the instrument or accessories to the manufacturer.

For decontamination of the RIDA® ABSORBANCE 96, follow the government guidelines for inactivation of organisms used in biological laboratories.

6.4 Repairs

Repairs on the device may only be carried out by the manufacturer. Please contact the service department. The product liability expires if the device is modified by unauthorized persons and by the installation of parts of a different type.

7 Performance data and technical data

Tab. 5: Performance data and technical data

Parameter	Value
Type of product	Absorbance reader
Product name	RIDA® ABSORBANCE 96
Service life	10 years at an average use of 4 h/day
Housing material	Aluminum
Place for use	Laboratory
Degree of contamination	2
Temperature (storage/measurement)	5 - 40 °C
Temperature (transport)	-10 - 50 °C
Relative humidity (storage/measurement/transport)	0 - 90 %
Air pressure (measurement)	900 - 1070 hPa
Air pressure (storage/transport)	600 - 1070 hPa
Height	Max. 5100 m
Measurement method	Absorbance
Light source	LED
Measuring technique	Endpoint and kinetic
Measurement channels	96
Filter	4
Wavelength	540, 450, 405, 630 nm
Measuring range	0 - 3.5 OD
Resolution	0.001 OD
Accuracy*	405 nm: 0.000 - 2.000 OD ≤ (1.5 % + 0.010 OD) 2.000 - 3.000 OD ≤ (3.0 % + 0.010 OD) ≥ 450 nm: 0.000 - 2.000 OD ≤ (1.0 % + 0.010 OD) 2.000 - 3.000 OD ≤ (1.5 % + 0.010 OD)
Reproducibility**	0.000 - 2.000 OD ≤ (0.5 % + 0.005 OD) 2.000 - 3.000 OD ≤ (1.0 % + 0.010 OD)

Linearity***	405 nm: 0.000 - 2.000 OD ≤ 1.5 % 2.000 - 3.000 OD ≤ 3.0 % ≥ 450 nm: 0.000 - 2.000 OD ≤ 1.0 % 2.000 - 3.000 OD ≤ 1.5 %
Connection to computer	USB 2/USB 3 with 5 VDC and max. 3 A
Dimensions	55 x 96 x 154 mm
Power supply	5 VDC
Nominal value/characteristic of the fuse	1 A / very fast-acting
Power input	2.5 W
Weight	900 g

*Accuracy is the maximum deviation between the determined value and the true value.

**Reproducibility is the maximum deviation between the determined values when the measurement is repeated directly.

***Linearity is the maximum deviation between the true and the determined increase of the value.

8 Warranty

The warranty period of the RIDA®ABSORBANCE 96 is 12 months. If a defect occurs in your device during this period, please contact our service staff directly. The RIDA®ABSORBANCE 96 may only be operated in technically perfect condition.

In the event of defects that could endanger employees or third parties, the RIDA®ABSORBANCE 96 may only be used again after it has been repaired. This warranty does not cover damage caused by improper use or external mechanical influences, transport damage or unauthorized intervention in the device by unauthorized persons.

Explanation of symbols (device nameplate)

General symbols of device nameplate:

Follow the instructions for use

Article number

Serial number

Manufacturer + Adress

CE mark

Disposal note/WEEE

Data matrix code

EC Declaration of conformity/CE mark

The RIDA®ABSORBANCE 96 is ROHS compliant according to Directive 2011/65/EU.

The RIDA®ABSORBANCE 96 was developed and tested according to CISPR 11 Class A. This declaration loses its validity, if the product is modified without our approval or if it is used for purposes other than those for which it was intended.

R-Biopharm AG

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64297 Darmstadt, Germany

Phone: +49 (0) 61 51 - 8102-0

Fax: +49 (0) 61 51 - 8102-40

E-mail: info@r-biopharm.de

www.r-biopharm.com

R-Biopharm AG

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64297 Darmstadt, Germany

Phone: +49 (0) 61 51 - 8102-0

Fax: +49 (0) 61 51 - 8102-40

E-mail: info@r-biopharm.de

www.r-biopharm.com

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© Copyright 2019 by R-Biopharm AG