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wild yeasts 1

1. Intended use

Detection of wild yeasts in beverages by real-time PCR.

The following yeasts are detected:

<i>Dekkera anomala</i>	<i>Dekkera bruxellensis</i>	<i>Dekkera custersiana</i>	<i>Dekkera naardenensis</i>	<i>Debaromyces hansenii</i>
<i>Hanseniaspora guillermondii</i>	<i>Hanseniaspora osmophila</i>	<i>Hanseniaspora uvarum</i>	<i>Issotchenkia orientalis</i>	<i>Kazachstania exigua</i>
<i>Kluyveromyces marxianus</i>	<i>Metschnikowia pulcherrina</i>	<i>Pichia anomala</i>	<i>Pichia fermentans</i>	<i>Pichia membranaefaciens</i>
<i>Saccharomyces cerevisiae</i> <i>var. diastaticus</i>	<i>Saccharomyces ludwigii</i>	<i>Torulasporea delbrückii</i>		

2. Test principle

The TaqMan® real-time PCR is based on hot-start-PCR and sequence-specific dual labelled probes (FAM/DQ) which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument.

To avoid false negative PCR-results an Inhibition Control (HEX/DQ) is amplified together in one reaction vessel with the specific sequence. The system contains dUTP. Optional: the use of Uracil-N-Glycosylase will eliminate any contamination with Uracil containing amplicons from former PCRs (the enzyme is not part of this kit).

3. Kit contents

The kit contains sufficient reagents for 50 reactions:

1 x	Premix	white cap
1 x	Dye Mix (freeze-dried, incl. IC-DNA)	dark vial, red cap
1 x	ddH ₂ O	colourless cap
1 x	Control-DNA (freeze-dried)	yellow cap

4. Storage conditions

The Dye Mix and the Control-DNA are freeze-dried, they have to be solved in ddH₂O prior to use (see 6.1).

Do **not** freeze the lyophilized Dye Mix and lyophilized Control-DNA.

The PCR reagents should be stored at 2 – 8 °C (35 – 46 °F).

Keep Premix for storage at - 20 °C (- 4 °F) after opening. Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The Dye Mix contains the fluorescent labelled probes and should be handled light protected.

All reagents are stable for 12 months, if they are stored correctly.

5. Materials required, but not provided

5.1. Instruments

Real-time PCR instrument with FAM- and HEX/JOE-channel

Centrifuge for 1.5 – 2.0 mL reaction tubes

Centrifuge for multiwell plates or stripes

Pipettes

“Vortex”

5.2. Reagents and plastic ware

sterile ddH₂O

sterile reaction vessels 1.5 – 2.0 mL

sterile optical tubes, plates, capillaries

sterile filter tips

optional: Uracil N-Glycosylase (0.01 U/μL added to the PCR reaction mix)

LC480: Use a Colour compensation prior to first use of the kit

6. PCR

6.1. PCR-Setup

When using the kit for the first time, the freeze-dried kit components have to be shortly centrifuged and carefully resolved:

- add 80 μL sterile ddH₂O to freeze-dried Dye Mix
- add 55 μL sterile ddH₂O to the freeze-dried Control-DNA
- after 15 minutes mix well

Before every use thoroughly mix all PCR components and centrifuge briefly.

PCR-reaction Setup:

PCR-components	amount (μL)
Premix	13.5
Dye Mix	1.5
Sample-DNA	2.5-5.0*
Total volume	20.0

* if using the Simplex Easy DNA- or QuickGEN Yeast Sample Preparation kit: add 2.5 μL DNA and 2.5 μL sterile ddH₂O

1. Prepare a mastermix by mixing Premix and Dye Mix.
2. Multiply said volumes with the number of PCR preparations including controls (positive control, negative control, extraction control), taking into account pipette reserves of approximately 5-10 %.
3. Divide 15 μL of the PCR-mastermix among the individual reaction vessels, making sure that, prior to the first filling, the tip of the pipette has been moistened.
4. Add 5 μL sample DNA, add 5 μL of the Control-DNA for the PCR positive control, add 5 μL of the extraction control and 5 μL of ddH₂O for the negative control* reaction. Use a fresh tip with each DNA filling.
5. Close the tubes immediately and centrifuge them shortly.
6. Place the tubes in the PCR-machine and start run.

Very important: * Please fill up the negative control with 5 μL ddH₂O to avoid unspecific amplification.

Work swiftly to avoid warming up and keep away from light.

6.2 PCR-Program

6.2.1 PCR-Program LC480

1. Click the button **tool** on the right side in the window **LightCycler 480 Software release 1.5.0. SP1**.
2. Click the button **Detection formats** at the left side of the menu bar.
3. Click **New** in the window **Detection formats** and name the experiment
4. Open the window **Filter Combination Selection** and choose the following filter combinations: 465-510 / 533-580.
5. Open the window **Selected Filter Combination** list and add the following amounts.

Excitation Filter	Emission Filter	Name	Melt factor	Quant factor	Max. Integ. Time
465	510	465-510 Fam	1	10	2
533	580	533-580 Hex	1	10	2

6. **Close** the window.
7. Click the button **New Experiment** on the right side of the menu bar.
8. From the pull-down menu **Detection formats** choose the defined experiment, click the button **Customize** and check the detection formats. All of them have to be activated.
9. Click the button **ok**
10. Define the following program:

1. Programm Name: **Heat**

Cycles **1** Analysis Mode **None**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)
37	None	00:15:00	4.40		0	0	0
95	None	00:15:00	4.40		0	0	0

2. Programm Name: **Ampli**

Cycles **40** Analysis Mode **Quantification**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)
95	None	00:00:15	4.40		0	0	0
62	Single	00:00:20	2.20		0	0	0

3. Programm Name : **Cool**

Cycles **1** Analysis Mode **None**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)
40	None	00:00:20	2.20		0	0	0

Optional: Saving the program as **run template**. The template button allows to select and apply a template to the currently open object and to save the currently open object as a template. Click the clamp beside the button **Apply Template**. Click **Save as template** and save the file in **templates**.

11. Click the button **Subset editor** on the left side
12. Click the button **+** and **New Subset 1** appears
13. Mark the wells in the **New Subset 1 Settings** window
14. Click the button **Apply**
15. Click the button **Sample editor** on the left side
16. **Very important:** Activate In the window **Step 1 Select Workflow: Abs.Quant**
17. Choose the subset **New Subset** in the window **Step 2 Select Samples**
18. Define your probes in the **Sample table**
19. Click the button **Experiment** and start the run with **Start run**

6.2.2 PCR-Program for other real-time machines

The Premix contains no ROX. This has to be done according to the settings of the real-time machine. **ABI 7500:** Under “Assign Targets and Samples” in ‘Select the dye to use as the passive reference’ choose „none“.

Mx3500P: Under Instruments in ‘Filter Set Gain Settings’ change the Filter set gain multiplier to CY5 x1, ROX x1, HEX-JOE x4 und FAM x8.

For the use of UNG the thermal cycler program has to be changed according to manufacturers` instructions.

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation	15 min	95 °C	
Cycling Denaturation	15 sec	95 °C	Cycle 40 x
Cycling Annealing/ Elongation	20-30 sec*	62 °C	

* please program 20 sec, if the real-time-PCR machine allows it

For MyGo Pro User :

Software 3.3.2

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation	15 min	95 °C	
Denaturation	15 sec	95 °C	Cycle 40 x
Annealing/ Elongation	20 sec	62 °C	

Advanced Settings

Integration time (s) 0.5
Acquisitions per Cycle 10

Software 3.4.8

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation	15 min	95 °C	
Denaturation	15 sec	95 °C	Cycle 40 x
Annealing/ Elongation	30sec	62 °C	

Advanced Settings

Integration time (s) 0.5

7. Evaluation

The evaluation has to be made according to the data analysis program recommended by the real-time instrument manufacturer.

For LC480: Activate Colour Compensation

wild yeasts 1-DNA: FAM-channel (LC480: 465-510 nm)


Inhibition Control-DNA: HEX-channel (LC480: 533-580 nm)

A sample is **wild yeasts 1 positive**, if there is a detectable fluorescence increase in the **FAM-channel (LC480: 465-510 nm)** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition-Control in the HEX-channel (LC480: 533-580 nm) may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). For negative controls it has to be positive.

A sample is negative, if there is no detectable fluorescence increase in the FAM-channel (LC480: 465-510 nm) and the positive controls have a positive fluorescence signal. The negative controls show no amplification. The Inhibition Control in the HEX-channel (LC480: 533-580 nm) has to be positive in the sample and in the negative controls, a false negative result due to inhibitory effects is then excluded.

analysis flowchart

LC480: analysis after activating Colour Compensation

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method. GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.