RIDASCREEN[®] SET A,B,C,D,E

Art. Nr.: R4101

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RIDASCREEN® SET A,B,C,D,E

Brief information

RIDASCREEN® SET A,B,C,D,E (Art. No.: R4101) is an enzyme immunoassay for the identification of staphylococcal enterotoxins A, B, C, D and E in fluid and solid foods as well as in bacterial cultures.

All reagents required for the enzyme immunoassay are contained in the test kit.

The test kit is sufficient for 12 determinations.

A microtiter plate spectrophotometer is required for determination.

Sample preparation: extraction according to a simplified method for working up of samples (partly applicable, s. point 9. „Sample preparation“) by homogenization with buffer and sample preparation

extraction according to the „Official European Preparation Method“ (dialysis concentration method, (see point 9.5.)

Time requirement: sample preparation for 10 samplesca. 1 h (simplified working up)

sample preparation for 10 samplesca. 19 h (official method, dialysis overnight)

test implementation (incubation time).....2 h 45 min

Limits of detection:

Simplified working up liquid samples 0.25 ng/ml Toxin
solid samples0.375 ng/g Toxin
supernatants from bacterial cultures 0.25 ng/ml Toxin

Dialysis concentration liquid samples 0.05 ng/ml Toxin
solid samples0.05 ng/g Toxin

The specificity of the RIDASCREEN® SET A,B,C,D,E test was determined by analyzing the cross-reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system

due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website www.r-biopharm.com/products/food-feed-analysis.

Related products

RIDASCREEN® SET Total (R4105; 96 determinations)

RIDASCREEN® SET Total (R4106; 48 determinations)

1. Intended use

RIDASCREEN® SET A,B,C,D,E is a sandwich enzyme immunoassay for the identification of staphylococcal enterotoxins (SET) A, B, C, D, E in fluid and solid foods as well as in bacterial cultures. Based on its sensitivity the RIDASCREEN® SET A,B,C,D,E test is consequently clear superior to the immunodiffusion procedure which has a detection limit of 0.1 mg/ml.

2. General


Staphylococci belong to the family of micrococci. *Staphylococcus aureus* and *Staphylococcus hyicus* produce one or more heat stable proteins which behave as enterotoxins. These are the causative agents for a number of food poisoning cases.

The main causative agent of food poisoning are enterotoxins of *Staphylococcus aureus* next to the Salmonella. Generally, it is assumed that a population of 5×10^5 cells of enterotoxin-producing *Staphylococcus aureus* strains per gram of food is required to lead to an intoxication. However, other studies show that only 100 - 200 ng of staphylococcal enterotoxins can lead to symptoms of food poisoning.

SET intoxications have been frequently associated with pasta, finished meat products, ham, pies, chicken meat products, fish, fish products, milk, milk products, ice-cream, egg products, salads, pastries and cake stuffing as well as preparations from these food products. The enterotoxins of the serological group A, B, C, D, E are very significant.

3. Test principle

The basis of the test is the antigen-antibody reaction. The wells A - E and H in the microtiter strips are coated with specific antibodies against staphylococcal enterotoxins A, B, C, D and E, the wells F and G serve as controls and are coated with antibodies of non-immunized animals in the following pattern:

Micro titer strips of the RIDASCREEN® SET A,B,C,D,E enzyme immunoassay

This sequence is obtained when the microtiter strips are inserted properly into the frame (narrow end of the strip → at the top, wide end of the strip → at the bottom). By adding the sample solution to the wells A to G (see the letters on the left side of the frame) and the positive control to well H, present toxins will bind to specific capture antibodies. Sample components not bound by the antibodies are then removed in a washing step. The immobilized toxins were bound by a mixture of specific antibodies conjugated to biotin. Any unbound conjugate is then removed in a washing step. After addition of Streptavidin-POD as well as an additional washing step the detection is performed by adding the substrate/chromogen solution. Bound Streptavidin-POD conjugate converts the reddish chromogen into a blue product. The addition of the stop solution leads to a color change from blue

to yellow. The measurement is made photometrically at 450/620 \pm 10 nm; the absorbance is proportional to the SET concentration in the sample.

4. Reagents provided

Each kit contains all reagents which are sufficient for 12 tests (from the 96m wells of the microtiter plate one strip of 8 single wells is necessary per single test).

Component	Cap color	Format		Volume
Microtiter plate	-	ready-to-use		96 wells
Positive control	red	ready-to-use	red stained	2 ml
Wash buffer	brown	Concentrate	10 x	100 ml
Conjugate 1	red	ready-to-use	green stained	11 ml
Conjugate 2	black	ready-to-use	blue stained	11 ml
Substrate/Chromogen Red Chromogen Pro	brown	ready-to-use	red stained	13 ml
Stop solution	yellow	ready-to-use		14 ml

5. Materials required but not provided

5.1. Equipment:

- analytical balance and weighing vessels
- Blender, Mixer or equivalent (for sample homogenization). For whey powder and all type of food difficult to mix, a turrax is highly recommended in order to have a homogeneous sample.
- 50 ml tubes
- optional: sterile filter cartridges
- 100 μ l micro pipette
- Incubator 35 - 37 °C (95 – 98.6 °F)
- Multi channel pipette or microwell plate washer
- optional: microwell plate spectrophotometer (450/620 \pm 10 nm)

Additional equipment for sample preparation according to the European official screening method (EOSM) Version 5, September 2010:

General remark: it is strictly recommended to use only laboratory vessels of glass-ware or polypropylene-ware (e.g. tubes, funnels, beakers, vials), to avoid adsorption of toxins

- Shaker for beakers (room temperature)
- Centrifuge, 3130 to 10 000 g, capable of being refrigerated at 4°C, centrifuge tubes
- Dialysis membrane, MWCO: 6 – 8 kD, flat width: 23 ± 2 mm (e.g. Spectra/Por[®]1, ref: 132 650, Spectrum)
- Closures, sealing width = 35 mm (e.g. Spectra/Por[®], ref: 132 736, Spectrum)
- pH-meter
- 50 ml tubes
- Funnel
- Glass-wool
- Tray
- Fridge (5 ± 3 °C/ 41 ± 5.4 °F) and freezer (≤ -18 °C/ ≤ -0.4 °F)
- Vortex
- Microwell plate spectrophotometer ($450/620 \pm 10$ nm)

5.2. Reagents:

- PBS buffer, pH 7.4, (0.55 g $\text{NaH}_2\text{PO}_4 \times \text{H}_2\text{O}$ + 2.85 g $\text{Na}_2\text{HPO}_4 \times 2 \text{H}_2\text{O}$ + 8.7 g NaCl, fill up to 1000 ml with distilled water)
- Distilled or osmosed water (optional: sterile water)
- n-heptane (for samples with high fat content)
- Brain-Heart-Infusion (BHI) for the pre-enrichment of potentially toxin forming Staphylococci strains. Ask your local producer/distributor of culture media for the supply with BHI

Additional reagents for EOSM:

- Hydrochloric acid (5N and 1N)
- Sodium hydroxide (5N and 1N)
- Polyethylene glycol 20 000 (PEG), quality for synthesis

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). **Do not freeze any test kit components.**

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label

The test kit can be regularly used at least up to the expiry date (indicated on the kit package), if stored correctly.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the reddish substrate/chromogen solution prior to test implementation
- a value of less than 1.0 absorbance units ($A_{450/620 \pm 10 \text{ nm}} < 1.0$) for the positive control

9. Preparation of Samples

The samples have to be stored at 5 ± 3 °C (41 ± 5.4 °F) until extraction. The samples should be completely defrosted before the extraction step.

Important notes: For good precipitation of solids and for complete separation of the phases while centrifugation it is necessary in some cases to increase centrifugation speed and/or time.

Food proteins may interact with toxins or antibodies of the test kit. In case of interaction they will affect the detection of the toxins tremendously. Therefore it is strongly recommended to prepare food samples according to the „Official European Screening Method of the European Reference Laboratory for Coagulase-positive Staphylococci” (see 9.5.).

The simplified methods for sample working up described under 9.1. to 9.3. are possibly not suitable for extraction of SET from difficult matrices (e.g. fish, chocolate, acidified vegetables/pickles). For these foods particularly low recovery rates as well as unspecific bindings of matrix proteins to the test antibodies have been observed.

A list of recovery rates for SET according to a variety of different food matrices is available on request.

9.1. Simplified working up of milk

- centrifuge milk samples (10 - 25 ml), especially raw milk samples, in cool condition: 10 min / 3500 g / 10 °C (50 °F)
(if no cooled centrifuge is available pre-cooling of samples is necessary)
- remove upper cream layer by absorbing it generously
- use 100 µl per well in the assay

9.2. Simplified working up of pasta and rice (cooked), meat, ice cream, processed food and other food with a fat content of less than 40 %

- mince 10 - 25 g of the sample and homogenize with 1.5 ml of PBS buffer (pH 7.4) per g sample (e.g. 10 g sample + 15 ml buffer)
- shake for 15 min
- centrifuge: 10 min / 3500 g / 10 °C (50 °F)
- if necessary, remove upper fat layer
- use 100 µl per well in the assay

9.3. Simplified working up of food with a fat content of more than 40 %

- mince 10 - 25 g of the sample and homogenize with 1.5 ml of PBS buffer (pH 7.4) per g sample (e.g. 10 g sample + 15 ml buffer)
- shake for 15 min
- centrifuge: 10 min / 3500 g / 10 °C (50 °F)
- transfer supernatant to another centrifugal vial, add the same volume of n-heptane and mix thoroughly for 5 min
- centrifuge: 5 min / 3500 g / 10 °C (50 °F)
- remove upper heptane layer generously, avoid to transfer heptane residues into the wells
- use 100 µl of the resulting aqueous (lower) phase per well in the assay

9.4. Bacterial cultures

Strains of the specie *S. aureus* (as well as *S. hyicus* or *S. intermedius*) which are potentially toxin forming have to be pre-enriched in Brain-Heart-Infusion (BHI) to ensure optimal formation of enterotoxins.

Important note: Prior to analysis please make sure that all strains to be tested are present in pure cultures

- centrifuge supernatants of microbiological fluid cultures 5 min / at a minimum of 3500 g / 10 °C (50 °F)
- sterile filtration of the supernatant is strongly recommended as not precipitated or re-suspended cells may influence the test reaction
- use 100 µl of the filtrate per well in the assay

Remark: Background effects which may occur (too high values for the negative controls) can be avoided, if the sterile filtrated supernatant is further diluted with PBS buffer.

This is also valid if OD-values measured for detected toxins are out of range of linearity (about ≥ 3.0 , please consider manufactures' instructions).

Sample preparations as well as culture supernatants can be stored at 5 ± 3 °C (41 ± 5.4 °F) if the analysis can be performed within 48 h. For longer storage periods the solutions have to be frozen at ≤ -18 °C (≤ -0.4 °F).

Frozen sample preparations or culture supernatants have to be completely defrosted and brought to room temperature before performance of analysis. Implement the once thawed samples into RIDASCREEN® SET A, B, C, D, E immediately. They can't be stored deep frozen again.

Sample preparations or culture supernatants prepared for analysis with RIDASCREEN® SET Total can be used directly for toxin identification with RIDASCREEN® SET A,B,C,D,E.

9.5. Sample preparation according to the European official screening method, Version 5, September 2010

9.5.1. Preparatory steps before toxin extraction

As staphylococcal enterotoxins can be heterogeneously dispatched in the sample, mix the whole sample if possible, or a representative part of it, with a mixer.

Weight 25 g \pm 0.1 g of the mixed sample and transfer this test portion into a beaker.

NOTE 1: In case of cheese with rind, take about 10% of rind for 90% of cheese.

NOTE 2: In case of powder samples, reconstitute the sample by weighting 12.5 g of sample and 12.5 g of distilled water or follow the manufacturer's instructions (ex: milk powder at 10%).

NOTE 3: In the case of a suspected Staphylococcal Food Poisoning Outbreak (SPFO), the minimal size of the test portion to perform the analysis is equal to 12.5 g.

9.5.2. Extraction step

Add 40 mL of warm distilled or osmosed water (38 °C \pm 2 °C/100.4 \pm 7.2 °F) to the test portion and homogenize the mixture by using a turrax, a blender or a stomacher. Rinse the system with distilled water.

NOTE 1: In case of liquid product, don't add 40 mL of distilled water.

Allow the toxin to diffuse by shaking the sample, at room temperature for at least 30 min.

Acidification step: Acidify the mixture with a few drops of hydrochloric acid in order to obtain a **pH between 3.5 and 4.0**.

*NOTE 2: During the acidification of the sample in order to keep enterotoxins in a good shape, a specific attention has to be given to: i) a pH-meter must be used and ii) **pH between 3.5 and 4.0** before centrifugation has to be respected.*

Be careful not to have a pH < 3.0 using hydrochloric acid. If the pH gets < 3.0, take another 25 g test portion and proceed as described on 9.5.1.

Centrifuge the mixture at least at 3130 x g for 15 min at **4 °C (32.9 °F) or room temperature**. Transfer the supernatant in a beaker.

NOTE 3: Do not hesitate to rinse with distilled water at each step to recover a maximum of toxin.

NOTE 4: If the supernatant is not clear enough, centrifuge again as described above.

*NOTE 5: The pH of the supernatant after the first centrifugation has to be < 4.5. If it is not the case, acidify until obtaining a **pH between 3.5 and 4.0** and centrifuge again as described above.*

Neutralization step: Neutralize the mixture using NaOH solution in order to obtain a **pH between 7.4 and 7.6**. Centrifuge again as described above. Recover the whole neutralized aqueous phase.

*NOTE 5: During the neutralization of the sample in order to keep enterotoxins in a good shape, a specific attention has to be given to: i) a pH-meter must be used and ii) **pH between 7.4 and 7.6** after neutralization has to be respected.*

Be careful not to rise above pH 9.0. If the pH is > 9.0, take another 25 g test portion and proceed as described on 9.5.1.

9.5.3. Dialysis-concentration step

For each sample:

- Prepare a 30% (w/v) PEG solution (30 g PEG / 100 mL distilled water).
- Cut about 50 to 60 cm of a dialysis membrane.
- Soak the membrane into distilled water following the instructions of the manufacturer (at room temperature during at least 30 minutes).
- Rinse the membrane with distilled water (outside and inside).
- Lock one end of the membrane with a closure, fill it up with the neutralised aqueous phase as prepared on 9.5.2 by using a funnel and a small piece of glass-wool to discard suspended particles. Lock the other end of the membrane with a second closure.

NOTE 1: If the sample to analyze contains high amounts of salt or sugar, conduct a dialysis under agitation with 2 L of distilled water, two times during one hour.

- Lay down the filled dialysis membrane in a tray containing the 30% (w/v) PEG solution.
- Allow the extracts to concentrate overnight at 5 ± 3 °C (41 ± 5.4 °F).

NOTE 2: If the extract is not concentrated enough, lay down it in the PEG solution for more time or add some powder of PEG.

- Take the dialysis membrane out of the PEG solution and rinse the outer-part of the membrane with distilled water to remove all traces of PEG.

Recover the concentrated extract using:

- the PBS solution in the case of extract containing milk or milk products**
- osmosis water in the case of extract containing no milk nor milk products**

Rinse well the inner-part of the dialysis membrane to obtain a final concentrated extract mass ranging from 5.0 g to 5.5 g (maximum 5.8 g for the stick extracts).

During this step, it is recommended:

- To rub the inner-parts of the dialysis membrane (one against another inner-part) with the fingers in order to take off and to recover the maximum of enterotoxins.
- To pour small drops of PBS or osmosis water (several additions) and to rinse in several times the inside membrane in order to recover all enterotoxins.

Transfer carefully the concentrated extract into a glass vial.

NOTE 3: If the initial weight of the sample to be tested is lower than 25 g (9.5.1, note 3, in case of SFPO), take care to obtain a final ratio equal to 5 between the weight of the concentrated extract and the weight of the test portion.

In case of SFPO or special the test portion mass may differ from the size of 25 g. The ratio of test portion mass to concentrated extract is shown in the following table:

Test portion mass	Concentrated extract
17.5 g – 25.0 g	3.5 g – 5.0 g (ratio 5:1)
12.5 g to 17.5 g	3.5 g (3.9 g max.)

NOTE 4: If the concentrated extract is analyzed within 48 h, store it at 5 ± 3 °C (41 ± 5.4 °F), otherwise store it at ≤ -18 °C (≤ -0.4 °F). The extract should be completely defrosted and homogenized before testing.

10. Test implementation

10.1. Preliminary comments

Bring all reagents as well as prepared samples to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **wash buffer** is provided as a 10fold concentrate. Before use, the buffer has to be diluted 1:10 (1+9) with distilled water (i.e. 10 ml buffer concentrate + 90 ml distilled water). Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (98.6 °F). The diluted buffer is stable at 2 - 8 °C (35 - 46 °F) for approx. four weeks.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

1. Insert as many microtiter strips into the microwell holder as samples are present. For each sample, all wells of a microtiter strip are required.
2. Add 100 µl of a sample to wells A to G of one microtiter strip, add 100 µl of the positive control to well H. Mix gently by shaking the plate manually. Cover the plate and incubate for 1 h at 35 - 37 °C (95 – 98.6 °F).
3. Dump the liquid out of the wells into a sink. Tap the microwell holder upside down onto a clean filter towel (three times in a row) to remove all remaining liquid from the wells. Fill the wells with 300 µl per well of wash buffer (10.1.) and pour out the liquid again. Empty the wells again and remove all remaining liquid. Repeat the washing step four more times.
4. Add 100 µl of conjugate 1 to each well. Mix gently by shaking the plate manually. Cover the plate and incubate for 1 h at 35 - 37 °C (95 – 98.6 °F).
5. Dump the liquid out of the wells into a sink. Tap the microwell holder upside down onto a clean filter towel (three times in a row) to remove all remaining liquid from the wells. Fill the wells with 300 µl per well of wash buffer (10.1.) and pour out the liquid again. Empty the wells again and remove all remaining liquid. Repeat the washing step four more times.
6. Add 100 µl of conjugate 2 to each well. Mix gently by shaking the plate manually. Cover the plate and incubate for 30 min at 35 – 37 °C (95 – 98.6 °F).
7. Dump the liquid out of the wells into a sink. Tap the microwell holder upside down onto a clean filter towel (three times in a row) to remove all remaining liquid from the wells. Fill the wells with 300 µl per well of wash buffer (10.1.) and pour out the liquid again. Empty the wells again and remove all remaining liquid. Repeat the washing step four more times.
8. Add 100 µl of substrate/chromogen to each well. Mix gently by shaking the plate manually. Cover the plate and incubate for 15 min at 35 – 37 °C (95 – 98.6 °F).
9. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450/620 ± 10 nm. Read within 30 minutes after addition of stop solution.

11. Results

For evaluation RIDA[®]SOFT Win.NET (Art. No. Z9996) which is a software especially developed for RIDASCREEN[®] enzyme immunoassays is available from R-Biopharm.

11.1. Quality control of the test (valid ranges for controls)

1. The absorbance value for the positive control should be equal or greater than 1.0.
2. The mean absorbance value for the negative controls should be equal or less than 0.2 units.

If these criteria have not been met, the following steps should be checked and corrected before repeating the test:

- check kit expiry date
- ensure sufficient time allowed for kit components to reach room temperature (20 – 25 °C / 68 - 77 °F)
- a new pipette tip must be used for each sample or control to avoid cross contamination
- check wash buffer for possible contaminations (poss. use sterile water for preparation)
- inadequate washing of wells/less washing steps as recommended (see 10.2.)
- contaminated pipettes (clean regularly)

11.2. Determination of threshold values

The mean absorbance value for the negative controls in wells F and G is determined separately for each sample. Add 0.15 to the mean value of the negative control to get the threshold value.

Example: Negative control 1 (well F) = 0.008
 Negative control 2 (well G) = 0.010
 Mean value = 0.009
 Threshold value = 0.009 + 0.15 = 0.159

11.3. Interpretation of results

1. A sample is considered POSITIVE for a certain SET if the test is valid (11.1.) and the absorbance measured at 450/620 ± 10 nm in the appropriate well is higher than or equal to the threshold value.
2. A sample is considered NEGATIVE for a certain SET if the test is valid (11.1.) and the absorbance measured at 450/620 ± 10 nm in the appropriate well is lower than the threshold value.

12. Crossreactivities of antibodies

Crossreactivities are known to occur between antibody/toxin: A/E, E/A, B/C and C/B. The percentage of crossreactivities is between 10 - 50%. Some food matrices (e. g.: milk, sausages) may enhance or weaken the crossreactivity. If crossreactivities are suspected the prepared samples should be diluted with PBS buffer (see below) and measured again.

To distinguish significantly between crossreactivity and actual detection the absorbance values of positive toxin measurements should be in a range between 1.0 to 2.5. For positive samples with higher extinction values than 2.5 the sample should be diluted, that the specific signal is in the range of 1.0 to 2.5.

Example 1 for crossreactivity between A/E:

A) Initial sample results:

Extinction value for detection of Toxin A: 4.1

Extinction value for detection of Toxin E: 3.6

B) The sample should be diluted to obtain an extinction value in the range between 1.0 and 2.5:

Extinction value for detection of Toxin A: 2.5

Extinction value for detection of Toxin E: **1.0**

Explanation: The signal of the sample in the cavity of toxin E is 1.0 and would be rated as strong positive. But it is only 40% of the specific signal of the sample in cavity A. This has to be evaluated as crossreactivity. The sample only contains toxin A.

Example 2 for no crossreactivity between A/E:

A) Initial sample results:

Extinktion value for detection of Toxin A: 4.1

Extinktion value for detection of Toxin E: 3.6

B) The sample should be diluted to obtain an extinktion value in the range between 1.0 and 2.5:

Extinktion value for detection of Toxin A: 2.5

Extinktion value for detection of Toxin E: **1.9**

The main difference between example 1 to example 2 is the extinktion value of the sample in cavity E (1.9), which is clearly above 50% of the specific signal of the sample in cavity A. In conclusion, there is no crossreactivity. The sample contains toxin A and toxin E.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

R-Biopharm AG

Postanschrift / Postal Address:

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64297 Darmstadt, Germany

Sitz / Corporate Seat: Pfungstadt

Tel.: +49 (0) 61 51 - 81 02-0

Fax: +49 (0) 61 51 - 81 02-40

E-mail: info@r-biopharm.de

www.r-biopharm.com

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are registered trademarks of R-Biopharm AG