Compact Dry AQ - medium for heterotrophic bacteria in water

100 plates/Platten/plaques/placas/lastre/placas

Art-No. HS9541

40 plates/Platten/plaques/placas/lastre/placas

Art-No. HS9542

English

Deutsch

Français

Compact Dry AQ is a ready to use, rehydrated plate with indicator for the detection and enumeration of heterotrophic bacteria in water

Specimen pretreatment**Viable count in water**

Drop 1 ml of the tested water on the middle of the Compact Dry plate.
The sample will diffuse equally all over the Compact Dry AQ plate.

Viable count using filter technology

For membrane water filter samples please use a 47 mm Mixed Cellulose Ester filter (MCE; 0.45 µm pore size).

Prior to applying the filter on the Compact Dry plate please pre-wet the plate with 1 ml buffer.
After water filtration place the filter in the center of the Compact Dry AQ Plate, grid side up. Roll top film down to minimize air bubbles or gaps between the filter and the Compact Dry AQ Plate.

Test instructions

1. Open the cap and drop 1 ml of specimen on the middle of the Compact Dry plate.
2. Specimen diffuses automatically and evenly into the sheet and transforms the dried sheet into a gel within seconds.
3. Put the cap again on the plate and write the information needed on the memorandum section.
4. Turn over the capped plate and put in the incubator in a horizontal position
5. After incubation count the number of all grown colonies underneath the plate. White paper placed under the plate helps to count the colonies.

Incubation time / Incubation temperature

36 °C ± 2 °C for 44 ± 4 h and
22 °C ± 2 °C for 68 ± 4 h (ISO 6222)
35 °C ± 1 °C for 48 ± 2 h (Filter/SMEWW Method)

Interpretation of the results

Colonies grown are almost all **red**. Yeasts tend to grow as **white-pink** colonies. Molds grow as typical 3-dimensional cottony form.

Storage and shelf life

Keep at room temperature (+ 1 to + 30 °C).
Total shelf life 24 months after manufacturing.

Notes

- High concentrations on plates (> 300 cfu) will cause the entire growth area to become reddish. In this case dilute the specimen.
- If the nature of sample does affect the reaction of the medium, inoculate the sample after the factor is eliminated by means of such as dilution and others. For instance; samples such as high viscosity, colored, reacting enzyme substrates, and too high or too low pH.
- After use please follow the current disposal regulations.
- The growth area is 20 cm². The back of the plate has a grid carved of 1 cm x 1 cm to make the colony counting easier. In case of any difficulties to count colonies due to large number of colonies grown, total viable count can be obtained by multiplying 20 by an average number of colonies per grid counted from several grids.
- Compact Dry plates are produced at an ISO 9001 certified site.

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